

## LAGIS EBH-0507 *FineRetriever™*

Features a push/pull rod and an introducer tube to facilitate bag opening and tissue capture during laparoscopic procedures.

In the fully deployed condition, the bag opening is maintained in a fully-open position by a metallic rim.

A string with a slip knot enables the closing of the bag after the specimen collection.



- 10mm diameter
- 300ml bag volume
- 5 in. x 7 in. bag size
- Strong, tear resistant
- Transparent pouch
- Radio-opaque rim



## EB-0507 *SimpleBag*

For quick, simple, safe specimen retrieval

PRODUCT CODE	DESCRIPTIONS
EBH-0507	<i>FineRetriever™</i>
EB-0507	<i>SimpleBag</i>

Packaging : 6 EA / Box